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## Editorial

Dear readers,

Turbulent times lie behind us and most likely ahead of us. The topic of hygiene has been on everyone's lips for months now. The most recent developments once again demonstrate that this topic has become an integral part of everyday life. All of us are having to confront it in our daily lives and even non-specialists now need to learn the correct procedures when it comes to taking protective measures and using disinfectants, for example.

This makes it all the more important that you, as experts, stay up to date with the latest developments and are informed about the topical issues surrounding hygiene and reprocessing. In his article "Mental hygiene in times of coronavirus – challenges and advantages", Dr Schimmelpfennig urges us to take care of our mental hygiene in order to get a better handle on crises and problems. Self-experiment is encouraged, especially at the present time. A critical look at our own behaviour in difficult situations and a certain degree of self-reflection can help us to find our inner balance.

In the hygiene and hospitals column, Ines Konschake shines a light on the topic of manual vs. machine-based reprocessing of TEE probes – a medical device that many hospitals only have in small quantities due to its high acquisition costs and that needs to be made available again very quickly due to the high number of patients requiring it: a real reprocessing challenge.

I hope you enjoy reading this issue of aseptica. Stay healthy!

Stella Nehr-Werner

[www.aseptica.com](http://www.aseptica.com)  
Extensive archive

## Report

### Rats wearing masks – artist Banksy spray paints the underground

Anonymous street artist Banksy has spray painted images of rats wearing masks on the interior and windows of trains on the London underground. "If you don't mask – you don't get." commented the artist on Instagram, posting a video of his work. The video shows what is likely to be the man himself, dressed in a protective suit, gloves, high visibility vest and, of course, a mask. In May, the artist donated one of his pieces to a hospital in Southampton. He left behind a note stating: "Thanks for all you're doing. I hope this brightens the place up a bit." The piece depicts a boy playing with a nurse doll. She looks as though she is flying in his hand, like a superhero.

The rats in the underground have since been removed by the London transport authority due to strict anti-graffiti regulations. The video has been watched more than four million times – so the message continues to be heard and shared.

Source: aerzteblatt

## Contents

### Latest News

Basic hygiene during the COVID-19 pandemic in dental practices 25

Mental hygiene in times of coronavirus – challenges and advantages 29

### Hospitals & Hygiene

*Pseudomonas aeruginosa*: a major threat to hospital water 32

### Info from Industry

ProCare Protein Test: the faster way of testing for proteins 34

New SL 2002 validation set for washer-disinfectors 34

### Technology & Hygiene

Manual or machine-based reprocessing of TEE probes: a comparison of advantages and disadvantages 35

Reprocessing thermolabile endoscopes – the latest information and discussions 44

Legal Notice 47

# Basic hygiene during the COVID-19 pandemic in dental practices

*Stella Nehr-Werner, Ulrike Weber*

Recent months have been consumed by the COVID-19 pandemic and have also posed huge challenges for the licenced dental sector. There has been a serious lack of personal protective equipment, uncertainty for staff and patients and a continuous stream of new information about the SARS-CoV-2 virus. The struggles caused by the virus and the need to process this new information are still present today. One thing has remained crucial at all times during this crisis and has provided a sense of safety as new measures are established – the principles of basic hygiene. Hygiene measures are based on the assumption that potentially infectious patients may enter the practice or be treated in the practice at any time.

The following will be an examination of the individual aspects of basic hygiene in the dental sector and conclusions will be drawn about the special circumstances of the SARS-CoV-2 pandemic. Specific recommendations for hygiene measures within practices are based on national recommendations, risk assessments and hazard evaluations and are not included in this article due to the highly dynamic nature of the situation and the speed at which information is changing. In Germany, the health authorities as well as the German Dental Association (BZÄK) and the different regional dental associations have specified special measures with regard to SARS-CoV-2.

## Basic hygiene in the face of coronavirus

As an enveloped virus, SARS-CoV-2 does not pose any major challenges in terms of hygiene measures when it comes to choosing products and processes. The efficacy of hand, surface and instrument disinfectants can also be determined quickly thanks to the principle of efficacy ranges for disinfectants (e.g. limited virucidal activity). Instrument reprocessing focuses on thermal disinfection, which also covers SARS-CoV-2.

However, a prerequisite for the implementation of basic hygiene is that all relevant microorganisms, and not just SARS-CoV-2 (as an enveloped virus), are considered. Classic bacteria such as *S. aureus*, *S. pyogenes* and fungus such as *C. albicans* as well as other virus strains such as HBV, HIV and influenza continue to be a threat during the events of COVID-19 and need to be taken into account as target organisms. Basic hygiene requirements when selecting disinfectants and procedures are therefore just as important as ever: limited virucidal activity, virucidal for final disinfection, bactericidal, yeasticidal, A0 value of 3000 for reprocessing of invasive medical devices. Practice hygiene and principles must not be limited solely to coronaviruses. The requirement to cover the mouth and nose area in public spaces and the principle of social distancing are new compared with measures for dealing with other infectious agents and are highly important, but they are not specifically addressed in this article. The treatment of Covid-19 patients is also not considered in this article. In this respect, we refer you to national recommendations and advice from professional bodies such as the German Working Group for Hygiene in Dental Medicine (DAHZ).

## Different hygiene measures

The respective facility-specific measures for implementing hygiene measures are defined in the hygiene plan. Information and dental-specific templates for a general hygiene plan are available from BZÄK.<sup>1</sup>

## Personal hygiene

Personal protection starts with correct hand hygiene and also involves wearing face masks covering mouth and nose and wearing gloves. Particular attention needs to be paid to using the right personal protective equipment for different tasks. In the medical field, employees can expect to encounter a whole host of different hazards when performing their duties.

## Authors

Dr Ulrike Weber  
Scientific Application & Sales Support  
Customer Segments & Solutions  
Business Unit Miele Professional  
Carl-Miele-Straße 29  
33332 Gütersloh, Germany  
ulrike.weber@miele.com  
www.miele-professional.de

Stella Nehr-Werner  
Global Clinical Affairs Manager  
Sirona Dental Systems GmbH  
Fabrikstr. 31  
64625 Bensheim, Germany  
stella.nehr-werner@dentsplysirona.com  
www.dentsplysirona.com

**Fig. 1:** Dental handpiece with 4-nozzle spray without air supply.

**Fig. 2:** Dental handpiece with 4-nozzle spray with an air-water mixture.

Hazard evaluations carried out in collaboration with an occupational safety specialist and the company's medical officer provide information about this. Every workplace with its various operations and associated hazards is assessed and then appropriate guidelines are provided. For example, in the unclean zone of the reprocessing room, the following tasks need to be carried out during instrument preparation: disposal of single-use instruments as well as pointed and sharp objects, cleaning and disinfection of instruments and elimination of the associated risk posed by biological agents. This results in the following protective measures, for example: Firstly, wearing puncture-proof gloves, which are also ideally resistant to cleaning agents and disinfectants. Secondly, disposing of pointed and sharp objects in appropriate containers to avoid injury. Thirdly, wearing protective gowns to protect workwear against contamination with infectious material and wetting. Fourthly, wearing visors or goggles to protect mucous membranes and eyes during cleaning tasks. In addition to basic hygiene, the following measures should be implemented to protect staff and patients:

- Goggles/protective visors/shields, face masks covering mouth and nose, gloves and protective gowns should be worn for the entirety of treatment of the patient. The barrier function of protective clothing is only guaranteed when it is fitted correctly and when it is changed as often as needed.
  - During pauses in treatment, the recommended minimum safety distance between employees should be maintained.
  - Every patient should be asked if they have had any COVID-19 symptoms in the past two weeks by telephone before they attend their dental appointment and again when they enter the practice.
  - There should be no physical contact when greeting a patient.
  - Patients should be instructed to disinfect their hands when they enter the practice and before they leave.
  - Every employee should wear face masks covering mouth and nose inside the dental practice, even when talking to one another.
  - Reception areas can be protected by a liquid-tight partition.
  - The number of people waiting should be kept to a minimum to maintain an appropriate distance between patients.
  - Treatment schedules and appointments should be arranged so that distance rules can be adhered to in the reception and waiting area.
  - People attending with adult patients should wait outside the practice.
  - Patients should be instructed to touch as few surfaces as possible. This applies to door handles, for example.
  - Patients should wear face masks covering mouth and nose when entering and inside the practice.
  - Toys should be removed from the waiting area.
  - Team meetings should take place regularly. This ensures that measures and routines can be discussed, questions can be answered and any changes can be made.
- In the case of employees at increased risk of severe illness from COVID-19, an assessment should be made based on local infection rates of whether or not they are at increased risk from direct contact with patients.

## Water and aerosol

Aerosols are used in dental treatments, primarily in the form of sprays (e.g. cooling dental drills). It is essential to use cooling systems to ensure that the heat generated does not damage healthy tissue, e.g. during cavity preparation. To this end, there is the option of either using process water from the treatment units or connecting external cooling media (e.g. bottle system or sterile saline solution). From a hygiene standpoint, there are two major factors to consider here:

1. The water supplied to the treatment units needs to be of drinking water quality. Particular attention needs to be paid to the quality of the process water in the treatment units, especially after long periods of not being used (e.g. due to practice closures). The water must not pose any risk to the patient. Suitable measures also need to be put in place to ensure that no biofilm builds up in the water unit. This can be achieved by using continuous disinfection systems or by regularly disinfecting the water systems. Rinsing procedures need to be established as part of the practice's daily routine before the start of treatment and at the end of the working day to prevent stagnant water. Additional rinsing procedures between patients ensure that any contamination caused by reflux is flushed out and that the waterways are all supplied with fresh water.<sup>7</sup>

2. However, it is not just patients who need to be protected against the infection potential of water with appropriate measures. The cooling water also poses an infection risk to practice personnel due to the fine nebulisation of the water during treatment due to the supply of air and the mixing with saliva. Mouthwash with antiseptic properties (e.g. chlorhexidine, hydrogen peroxide) can reduce oral flora and therefore reduce pathogen exposure in sprays.<sup>2</sup> High-volume and efficient extraction systems can also reduce the dispersion of aerosols during treatment. Reducing the spray by decreasing or switching off the spray air minimises the risk thanks to reduced aerosol formation. It should also be mentioned that a dental turbine tends to produce more aerosol than an electrically driven high-speed motor due to the higher rotational speeds.

## Face masks covering mouth and nose

Person-to-person transmission via droplets (e.g. when talking, coughing, sneezing) is the main form of transmission

for coronaviruses, particularly for SARS-CoV-2. Due to the complex supply situation with regard to breathing masks, various methods for reprocessing breathing masks (primarily FFP2 and face masks covering mouth and nose) have been discussed. As these products are initially not developed and manufactured as reprocessable, not only do different reprocessing methods (e.g. damp heat, dry heat, UV, H<sub>2</sub>O<sub>2</sub>) need to be considered, but also the filtration, particle retention and fit properties after reprocessing. Different recommendations and information are being circulated in Europe and internationally. As a result, it can be concluded that reprocessing is the responsibility of the individual and must be subject to a very critical risk-benefit assessment. Due to the significantly increased availability of breathing masks (FFP2 and face masks), there is no longer any need for reprocessing.

## Clothing


The protective clothing used is either available as single-use items or multi-use items. With single-use clothing, it is important to adhere to correct disposal measures to ensure that the wearer does not come into contact with infectious or otherwise harmful material during disposal. If suitable personal protective equipment is reprocessed, observe the manufacturer's instructions, particularly with regard to use of cleaning agents and procedures as well as disinfectants and disinfection procedures. Workwear also needs to be subjected to a disinfecting reprocessing procedure.

## Surface hygiene

Surfaces may become contaminated through direct contact with hands or due to the descent of particles containing aerosols. The required frequency of disinfection procedures can be found in the cleaning and disinfection plan and the application times of the respective products must be observed. Suitable agents need to fulfil the bactericidal, yeasticidal and limited

Person-to-person transmission via droplets (e.g. when talking, coughing, sneezing) is the main form of transmission for coronaviruses, particularly for SARS-CoV-2.

**Fig. 3:** Instrument cycle.

virucidal activity spectra as a minimum requirement. Alongside the manufacturer information, the usual lists (Association for Applied Hygiene (VAH) and German Industrial Association for Hygiene and the Protection of Surfaces in Industrial and Institutional Applications (IHO)) also provide information on this. When selecting products, it is essential to consider material compatibility with various practice equipment surfaces and materials as well as the activity spectrum. Over recent months, chloride-containing products have also been frequently used due to supply shortages and the effectiveness, stability and material compatibility of these products were not always guaranteed.

### Reprocessing

The reprocessing of medical devices follows the standard, facility-specific instrument cycle principles (see Fig. 3). Machine-based reprocessing is the preference in this case.

### Summary

Additional coronavirus measures are being introduced based on the special transmission pathways for SARS-CoV-2 and a critical risk-benefit assessment. As a basic principle, hygiene measures need to be carefully attuned to the presence of germs just as they were before. Due

to the novelty of the virus and its special forms of transmission, recommendations are constantly being updated based on scientific findings. It is to be expected and hoped that the resulting measures will be more cautious and take into account risk and benefit assessments.

### References

1. German Dental Association (BZÄK) and German Working Group for Hygiene in Dental Medicine (DAHZ). HYGIENEPLAN / ARBEITSANWEISUNGEN für die Zahnmedizin (HYGIENE PLAN / WORK INSTRUCTION for dental medicine). Last edited 14.02.2020.
2. German Dental Association (BZÄK): Risikomanagement (Risk management), <https://www.bzaek.de/berufsausuebung/sars-cov-2/covid-19/risikomanagement.html>.
3. ECDC. COVID-19 infection prevention and control for primary care, including general practitioner practices, dental clinics and pharmacy settings 9 June 2020.
4. German Working Group for Hygiene in Dental Medicine: SARS-COV-2 / Covid-19 – Risikomanagement in Zahnarztpraxen (SARS-COV-2 / Covid-19 – Risk management in dental practices), <http://dahz.org/wp-content/uploads/2020/04/DAHZ-Stellungnahme-Corona-20.04.2020.pdf>.
5. RKI Richtlinie aus 2006 – Infektionsprävention in der Zahnheilkunde – Anforderungen an die Hygiene (RKI guideline from 2006 – Infection prevention in dentistry – Hygiene requirements).

# Mental hygiene in times of coronavirus – challenges and advantages

*Markus Schimmelpfennig*

When I first wrote about mental hygiene for a “specialist magazine for hospital and practice hygiene” (as this magazine defines itself) in the 01/2018 edition of *aseptica*, I was fully prepared for angry letters from readers – why on earth include such “psychological nonsense” in a scientific specialist magazine. But there was no such protest and some readers even expressed their approval, and this apparently encouraged the editors to ask me to make another contribution on this topic. And I was more than happy to oblige as I am still convinced that mental hygiene is a fundamental tool to help everyone with the practical implementation of hygiene.

Let's remind ourselves of the seven steps of mindfulness according to Kabat-Zinn:

1. Non-judging
2. Patience
3. Beginner's mind
4. Trust
5. Non-striving
6. Acceptance
7. Letting go

This is certainly a challenge in everyday life and requires constant practice, but it only works when you fully embrace these seven steps.

## Non-judging

Why should I be non-judging (except in an emergency, of course)? Because my immediate evaluation may cloud my view of the full picture and therefore lead me down the wrong path.

## Patience

Why should I be patient? Because grass does not grow quicker when I pull it and nothing can make a wound heal faster beyond good treatment and the patient's own self-healing powers, even if I want it to. Impatience rarely achieves anything, but drains a lot of energy.

## Beginner's mind

And why should I, as an “old hand”, have a beginner's mind when I am proud of being able to immediately assess a situation and take the necessary steps on the basis of my experience and knowledge, and rightfully so?

Because I can make critical mistakes if I am not prepared to at least occasionally question my own understanding and actions: “In this particular case, which looks the same as 1000 others that came before it, could I possibly be making the wrong assumption?”

## Trust

Why should I trust? In many respects, trust is the foundation of how I act and plan. If I feel as though I cannot rely on anything or anyone, not even on myself, I become incapable of action and lose all perspective. Of course, trust can be taken for granted and abused; most of us will have experienced this at least once. Nevertheless, trust is needed for us to have the courage to face life and be able to take action. If I want to feel solid ground beneath my feet, I need to be prepared to trust – maybe just a few individuals, maybe even just my dog to start with, but nothing can happen without trust.

## Author

Dr Markus Schimmelpfennig  
Specialist in public health, hospital hygiene,  
Marienkrankenhaus hospital Kassel  
Marburgerstraße 85  
34127 Kassel, Germany  
m.schimmelpfennig@marienkrankenhaus-kassel.de

Impatience rarely achieves anything, but drains a lot of energy.

“ But mental hygiene undoubtedly helps those who practise it to stay in touch with reality, to be better equipped to withstand any kind of stress, to maintain or regain a better balance, and to live a more balanced and contented life overall.

### Non-striving

Why should I be non-striving? Because striving requires a lot of energy and yet almost always results in failure. There are some things that I simply cannot strive for, such as love and healing. I can put everything into an endeavour, but if I don't succeed, I can preserve myself and my energy if I am able to accept my failure.

### Acceptance

And why should I accept and what should I accept? We should not, of course, accept everything: for example, unacceptable behaviour, injustice, obvious wrongdoing. Acceptance here means accepting something that I cannot change because it is out of my control without bearing any ill will. Feeling angry about the third red traffic light in a row when I am in a rush is pointless as I have no power to change the traffic light sequence – I may as well save the energy I would have expended during my temper tantrum for something more meaningful. What's more, acceptance is often the fundamental prerequisite for change. If I continue to see myself through rose-tinted glasses or deny that something is the way it is, there's no need for me to do anything about it. If I refuse to admit that I'm overweight and just make jokes about being “big boned”,

I'm not forced to question my lifestyle, and if I don't want to admit that I don't do enough exercise (which is not a problem for many of us because our work means that we are always on the go), I can continue living as a couch potato with a clear conscience. This attitude prevents acceptance becoming a driver for change.

### Letting go

And the final and perhaps most difficult step: Why should I let go and what should I let go of? I should let go of everything inevitable without inwardly battling with it. I find this extremely difficult, but if I don't let go, I will lose perspective and have no course of action. Grief is an example of this: lost loved ones will not come back to us no matter how much pain we feel. I have to face this pain, endure it, and then come to a point where I can let go of it, only keeping hold of the good memories, no matter how hard this might be. This is the only way I can feel solid ground beneath my feet again and start looking towards a new future.

And what does all this have to do with the coronavirus pandemic?

A great deal, as you will see. The coronavirus pandemic is an unintentional lesson in how much our mental hygiene can help us to cope in a bad situation, provided we are willing to practice and apply it.


Fig. 1: Seven steps to mindfulness.

Let's take a closer look: Right at the start, observant Chinese doctors identified a new acute respiratory disease caused by a novel virus, but they did not make judgements about the situation right away ("This can't be true", "If we report this, we'll get into a lot of trouble", "We have to keep quiet about this to protect our economy", etc.). They made their objective observation of a new acute respiratory disease caused by a new virus Sars-Cov-2 public to the world. The scientific community conducted extensive research worldwide and continues to do so in order to find out more about the virus and the disease it triggers and to develop and make available treatment and vaccination options as quickly as possible (= patience), but without trying to accelerate the process (= impatience = botched job = striving). This relies not only on tried-and-tested methods, but also on completely new approaches, e.g. in vaccine development (= beginner's mind). We then have to trust that everything humanly possible will be done to control this global threat, both scientifically and politically. However, we have also seen how a lack of acceptance worldwide ("it's not that bad really", "we won't put any unnecessary restrictions in place, we can handle it without impinging on freedom"), a lack of self-criticism ("I'm just great, everything I do is good!" as expressed by two heads of state on the American continent) and "rush jobs" = impatience ("Wouldn't the simplest thing be to just spray people with disinfectant?," also expressed by an American politician) can damage this trust.

We also need to take a good look at ourselves: The pandemic has been affecting us for about four months now with no end in sight and we are experiencing growing impatience within ourselves and in the people around us, a decreasing willingness to accept the necessary restrictions and a tendency to simply slip back into or even strive for the former freedoms we used to have without any scientific basis to back this up. But those who no longer perceive (without judgement) the danger, no longer comply with (accept) the necessary protective measures, whether that is in public or in private, and take back (strive for) their freedoms, will pay the price when infection rates rise again. And when confronting the pandemic, whether that is in terms of the world community or the fate of individuals, the limited amount of control we have and the associated feeling of powerlessness also (sadly) require us to let

go, whether that relates to the necessary economic restrictions or the patient who dies despite every effort or who could not be treated properly due to a lack of manpower or resources.

### Mental hygiene alone is not enough

This last point needs to be a trigger, both nationally and internationally, for a thorough analysis of the shortcomings in the health system with the aim of eliminating them as quickly as possible. This also involves addressing globalisation and equitable distribution problems (one world – one health) because in a time when a virus can travel from one corner of the world to every other corner of the world within 48 hours, we cannot be indifferent to who produces various medications and medical devices, where they are produced and under what conditions and how many of them are stockpiled and where because the profitability aspect alone is both a poor and incorrect benchmark for global and public health.

### Summary

The practical pursuit of mental hygiene brings benefits to anyone who wants to have a positive influence on their life and get a better handle on crises and problems. Mental hygiene cannot and is not intended to replace necessary measures in societies, states and the global community that go beyond the individual.

The coronavirus pandemic poses a significant, albeit involuntary, challenge for both the individual and the global community. But mental hygiene undoubtedly helps those who practise it to stay in touch with reality, to be better equipped to withstand any kind of stress, to maintain or regain a better balance, and to live a more balanced and contented life overall. Try it out – see for yourself how it can help you, especially during the coronavirus pandemic.

**Fig. 2:** The new normal, masks in the community.

# Pseudomonas aeruginosa: a major threat to hospital water

## Author

Dr Regina Pietsch  
Hospital hygienist  
Department of Hygiene and Infection  
Prevention, Mainz University Medical Center  
Hochhaus am Augustusplatz  
55131 Mainz, Germany  
regina.pietsch@unimedizin-mainz.de  
www.unimedizin-mainz.de/ahi

### Regina Pietsch

Section 3 of the German Drinking Water Ordinance<sup>1</sup> defines drinking water from various perspectives as a consumable that comes into contact with the human body. It is water used in the kitchen for preparing food and beverages, for clean-

ing dishes, for personal hygiene and for washing clothes and other laundry. This definition requires that, in accordance with Section 4 of the Ordinance, drinking water is free of pathogens that can adversely affect human health. *Pseudomonas aeruginosa* is one of these pathogens.

This germ belongs to the group of facultative human pathogens which can sometimes lead to severe illnesses not only as a result of drinking the water containing the germ, but also due to other contact with the water. This means *Pseudomonas aeruginosa* is a particular concern with regard to indoor piping systems, especially in large buildings and bathing water. This requires particular consideration in the medical field with the treatment of patients with pre-existing conditions. Outside of this field, the German Drinking Water Ordinance only requires proof of the germ if events make this necessary and does not specify a limit value for this. An exception to this is the rare occurrence that drinking water is poured into sealed containers. In this case – as specified in the Mineral and Table Water Ordinance<sup>2</sup> for bottling water from natural sources – the limit value of 0 germs in 250 ml applies. For bathing water, DIN 19643-3 stipulates the limit value of 0 germs in 100 ml.

### Germ with self-protection mechanisms

*Pseudomonas aeruginosa* is an unassuming germ that exists universally in soil and

water. Its reproductive capability spans over a wide temperature range of 4 °C – 44 °C (Fig. 1). This means that this germ can occur in both hot and cold water systems. It has the ability to form biofilms on surfaces. These protect the germ from running water and any disinfectants in the water. Expensive and laborious disinfection measures are often required to remove the germ from piping systems. Its defence mechanisms include a self-produced thick layer of mucus, which protects bacteria in the body against attack from antibiotics and the defences of the body's own immune system (Fig. 2). In addition to this, almost all *Pseudomonas aeruginosa* strains have intrinsic antibiotic resistance and can acquire additional resistance through gene transfer, making individual strains multi-resistant (3MRGN/4MRGN).

### Diverse symptoms

*Pseudomonas aeruginosa* can result in serious illnesses in hospitals, in retirement and care homes or for immunosuppressed patients who live at home.<sup>4</sup> In medical facilities, patients with catheter systems, receiving antibiotics or chemotherapy or with a chronic obstructive pulmonary disease (COPD) are particularly at risk. There is a wide range of illnesses that can occur, including urinary tract infections, wound infections, pneumonia and septicemia.

Being infected with *Pseudomonas aeruginosa* is particularly dangerous for patients with cystic fibrosis, in both a hospital setting and at home.<sup>5</sup> Outside of medical facilities, otitis externa, chronic wound infections, keratitis and folliculitis are typical infections caused by this germ, particularly via bathing water.<sup>4</sup>

### Infection risks in hospitals and private homes

The significance of evidence of *Pseudomonas aeruginosa* in the drinking water system of a building varies depending on its use. A distinction is therefore made between hospitals, retirement and care homes, hotels and residential buildings. It is important to note that basic mains water tests in accordance with the German

All of this proves that a hospital's water system is a complex matter that needs to be taken seriously due to the potential risk to patients.

Drinking Water Ordinance that detect *Escherichia coli* and coliform bacteria and that determine colony numbers provide no indication of contamination with *Pseudomonas aeruginosa*. This requires a specific test.

Numerous case reports have found that the germ can reach patients in hospitals via the piping system and result in infection.<sup>6</sup> This is particularly evident in intensive care units where various measures taken for the patient, such as body washing and oral hygiene, can lead to transmission of the germ. Predisposing factors for infection include the number of days spent on a ventilator for intubated patients and the duration of treatment with antibiotics that are ineffective against *Pseudomonas aeruginosa*. This germ can also result in infection in a domestic setting. Bacteriological monitoring of domestic plumbing systems is particularly important when a patient is immunosuppressed. Since the German Drinking Water Ordinance focuses on protecting the healthy population, the protection of patients with immunodeficiency is not sufficiently guaranteed without specific germ detection. Additional mains water tests need to be performed in order to prevent infection via water.

### Drinking water hygiene plan


The World Health Organization (WHO) has developed the concept of the water safety plan (WSP) to ensure a safe water supply.<sup>7</sup> This concept involves hospitals developing a drinking water hygiene plan in order to keep any bacterial contamination with *Pseudomonas aeruginosa* and other problematic germs under control. This requires the creation of a water hygiene room log and performance of a risk assessment. This plan will define sampling locations, the frequency of testing and measures to be taken in the event of suspicious lab results. A rinsing plan is required to prevent contamination of the water system, particularly in areas with low water usage. Regularly rinsing the pipelines prevents stagnation and germ growth in the pipelines and faucets. Technical maintenance of piping systems needs to be carried out by specialists in accordance with all recognised technical regulations.

All of this proves that a hospital's water system is a complex matter that needs to be taken seriously due to the potential risk to patients. Sterile filters on faucets and

shower heads in the event of positive findings can only be used as an acute measure. It does not replace identifying the cause and eliminating it, which is not an easy task with *Pseudomonas aeruginosa*.

**Fig. 1:** Differentiation of *Pseudomonas aeruginosa*: Fluorescence under UV light on cetrimide agar.

**Fig. 2:** Clear mucus formation around the *Pseudomonas aeruginosa* colonies beyond the filter.

### References

1. German Drinking Water Ordinance as published on 10 March 2016 (Federal Law Gazette (BGBl.) I p. 459), last amended by Article 1 of the Ordinance of 20 December 2019 (Federal Law Gazette (BGBl.) I p. 2934).
2. German Mineral and Table Water Ordinance of 1 August 1984 (Federal Law Gazette (BGBl.) I p. 1036), last amended by Article 25 of the Ordinance of 5 July 2017 (Federal Law Gazette (BGBl.) I p. 2272).
3. DIN 19643-1/2/3/4 Treatment of water of swimming pools and baths, publication date 2012-11.
4. Exner M, Nissing W, Behringer K et al.: Gesundheitliche Bedeutung, Prävention und Kontrolle Wasser-assoziiertes *Pseudomonas aeruginosa*-Infektionen (Significance to health, prevention and control of water-associated *Pseudomonas aeruginosa* infections). Hyg Med 2016; Suppl. 2: 3-27.
5. Von der Hardt H, Schwarz C, Ullrich G: Erwachsene mit Mukoviszidose – Es geht um mehr als die Lebensdauer (Adults with cystic fibrosis – It is about more than life expectancy). German Federal Health Bulletin 2012; 55: 558-567.
6. Mena KD, Gerba CP: Risk Assessment of *Pseudomonas Aeruginosa* in Water. Rev Environ Contam Toxicol 2009; 201: 71-115.
7. Schmoll O, Bethmann D, Sturm S, Schnabel B: Das Water-Safety-Plan-Konzept: Ein Handbuch für kleine Wasserversorgungen (The water safety plan concept: a manual for small water supplies). 2014. <https://www.umweltbundesamt.de/sites/default/files/medien/374/publikationen/wps-handbuch-web.pdf>.

## ProCare Protein Test: the faster way of testing for proteins

Fig. 1: Rapid protein test for routine monitoring.

Any proteins adhering to the load items are denatured and affixed to surfaces during a full cleaning and disinfection programme which includes thermal disinfection. At the end of a programme, an objective examination of cleaning performance on instruments by taking a sample and using analytical methods would no longer be possible. With its ProCare Protein Check, Miele Professional is presenting a rapid test which enables the simple and fast detection of protein residue on instrument surfaces visible to the human eye after thermal disinfection has taken place.

**Miele**

## New SL 2002 validation set for washer-disinfectors

Fig. 1: Set SL 2002.

The SL 2002 validation set from Xylem Analytics Germany Sales GmbH & Co. KG, ebro®, contains a temperature logger for measuring in the drain of the washer-disinfector and two temperature loggers with 4 sensors at different levels for measuring at the instruments to be cleaned. A pressure and temperature logger for rinsing pressure detection and temperature measuring and a logger for measuring conductivity are also included in the set. In combination with the read-out station and TÜV-certified software Winlog.validation, the set allows the user to start thermal testing for validation in accordance with ISO 15883 immediately.

**-ebro-**  
a xylem brand

# Manual or machine-based reprocessing of TEE probes: a comparison of advantages and disadvantages

*Ines Konschake*

Reprocessing probes for transoesophageal echocardiograms (TEE probes) and endoscopes is a major challenge in day-to-day work. Over the last few years, for example, there have been repeated reports of inadequate disinfection of medical devices like these. According to international reports, there have been instances in recent years of poorly reprocessed instruments being responsible for pathogen transmission.

TEE probes are classified as being in the “semi-critical B” risk group. This group is for medical devices that come into contact with mucous membrane or pathologically affected skin. There are more stringent requirements when it comes to reprocessing of TEE probes, as is the case for all medical devices in the “semi-critical B” category. Despite this, TEE probes that have not been properly disinfected are still used, which results in cross-contamination between patients. This is a serious health risk and a hazard for patients. A possible cause of this might be to the different reprocessing methods used.

## Requirements for reprocessing TEE probes

Two different methods for reprocessing TEE probes are available and common practice: manual reprocessing and automated reprocessing. Both methods are used in the day-to-day operation of hospitals whereas the manual TEE probe reprocessing is the most common. One reason for this is that there are currently only a few automated reprocessing methods available and the existing, well-established processes like manual reprocessing in hospitals cannot be changed overnight. In addition, TEE probes are also classified as medical devices with a risk level of “semi-critical B”. This means more stringent cleaning

and disinfection requirements for the reprocessing.<sup>1</sup> Despite this, day-to-day hospital operations indicate that shortcomings in reprocessing can occur due to insufficiently cleaned and disinfected TEE probes, therefore posing a contamination risk for patients.<sup>2</sup> In the case of manual reprocessing procedures, correct reprocessing always depends on the person completing the task. Due to the lack of reproducibility in the individual reprocessing steps by person-dependent application, errors cannot be excluded. However, valid efficacy requires proof of the constant effectiveness of the manual reprocessing method. According to the judgement of the Higher Administrative Court (OVG) of North Rhine-Westphalia, ref. 13 A 24222/09 of 29.09.2010, this means that manual reprocessing methods are not considered valid.<sup>3</sup> A recommendation from the German Society for Hospital Hygiene (DKGH) in April 2017 came to a similar conclusion.

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
Fig. 1: Probe disinfection plan.

**Fig. 2:** Risk classification of TEE probes.

“If multiple reprocessing options are specified, the safest method must be used. In general, this is the automated reprocessing in a washer-disinfector.”<sup>2</sup> The following product comparison shows three possible reprocessing methods for TEE probes and examines which method proves itself to be particularly effective, easy to perform and safe for patients and users.

Manual and automated TEE probe reprocessing using Tristel Trio Wipes, the Optima C10 and the Soluscope Serie TEE

One manual and two automated reprocessing methods were selected for the product comparison:

- Tristel Trio Wipes System (manual cleaning and disinfection)
- Optima C10 (semi-automated disinfection)
- Soluscope TEE (automated cleaning and disinfection)<sup>4</sup>

Next to the widely used manual immersion disinfection method, where the TEE probe is placed in a tube or tub with a disinfectant for a certain period of time, impregnated disinfection wipes, such as the Tristel Trio Wipes System, is another frequently used manual

reprocessing method. This manual wipe cleaning and disinfection method is often used to facilitate rapid reuse. The two other methods either involve a fully automatic and standardised reprocessing procedure, like with the Soluscope Serie TEE, or a semi-automated procedure for disinfection and rinsing, like with the Optima C10. All automated reprocessing are documented automatically on both devices.

The reprocessing procedure in detail: Tristel Trio Wipes, Optima C10 and Soluscope TEE

Pre-cleaning is a manual step in all three methods and takes place in the examination room after the protective cover has been removed from the TEE probe. The reprocessing procedure itself takes place in a separate reprocessing area, which has a clean zone and an unclean zone. The entire reprocessing procedure for each method is being compared, which comprises the following steps:

1. Cleaning
2. Disinfection
3. Rinsing
4. Drying

The duration of reprocessing cycle varies for the three methods. While the process duration using the Tristel Trio Wipes System (cleaning/disinfection/rinsing) always depends on the degree and extent of the contamination and the associated cleaning of the probe using the wipes and the short application time of the disinfectant (30 sec.), the semi-automated process (disinfection/rinsing) with the Optima C10 takes a total of 30 minutes and the fully automated Soluscope TEE (cleaning/disinfection/rinsing/drying) takes only 14 minutes. The Tristel Trio Wipes System is versatile and used for reprocessing various medical devices, such as transvaginal and transrectal ultrasound probes, tonometer prisms, etc. Optima C10 is also suitable for transvaginal and transrectal probes, while the Soluscope Serie TEE is exclusively for reprocessing TEE probes. A special feature of the Optima C10 is that an electrical residual current measurement (leak test) which is performed during reprocessing. The Soluscope Serie TEE also offers an automatic leak test with every reprocessing cycle. When using a manual method, the leak test needs to be performed prior to the actual reprocessing.

### Comparison of individual process steps for TEE probe reprocessing

#### Tristel Trio Wipes System

##### Scope of delivery:

The system comprises three wipes intended for the individual reprocessing steps:

1. Cleaning with the Tristel Pre-Clean Wipe
2. Disinfection with the Tristel Sporicidal Wipe and Activator Foam
3. Rinsing with the Tristel Rinse Wipe

##### Manual TEE probe reprocessing method:

The three reprocessing steps are standardised by the process sequence. The Tristel Sporicidal Wipe has a virucidal efficacy in accordance with EN and the Robert Koch Institute (RKI)/German Association to Combat Viral Diseases (DVV) and is listed in the Association for Applied Hygiene (VAH) list of disinfectants. There is also a detailed instruction for use available. A log book is included for documentation of each process step. The leak test is not part of the reprocessing procedure and needs to be performed separately.

#### Optima C10

##### Connection and operation:

The Optima C10 is connected to a standard power plug. The device is also connected to an existing water connection or drain. The disinfectant vapours need to be extracted via a direct connection of the Optima C10 ventilation filter to the room ventilation system or via an outside air connection. This minimises exposure to disinfection vapours. The device is also equipped with a filter system for reprocessing the water required for the rinsing process. The machine has a self-disinfection programme. The device is operated and various settings can be made via a colour touch screen display. Further information on a connectable barcode scanner for recording TEE probe identification and user data in the process documentation is not covered in the following presentation as this information is not available in either the manual<sup>5</sup> or the Optima C10 training video.

Fig. 3: TEE probe with protective cover.

*The semi-automated TEE probe reprocessing method:*

Automated reprocessing using the Optima C10 is limited to disinfection, the leak test (in accordance with VDE 0751) and rinsing and is documented via the integrated printer. The probe must be cleaned manually in a separate process before being placed in the device.

**Soluscope Serie TEE**

*Connection and operation:*

A standard power plug, a water connection and a drain are required for connection of the Soluscope Serie TEE. Connection to the room ventilation system is not required as no vapours are emitted during the process. The device is also equipped with a process-monitored HEPA filter and water filter for the water required during rinsing. A self-disinfection programme is another safety feature. The device is operated via a display. The device is certified as an automated reprocessing unit in accordance with ISO 15883-1.

*The automated TEE probe reprocessing method:*

The TEE probe reprocessing method of the Soluscope Serie TEE includes fully automatic performance of the three reprocessing steps: cleaning, disinfection and rinsing.<sup>6</sup> Drying also takes place inside the device. The leak test to identify damaged TEE probes is performed automatically with every reprocessing cycle. All process

parameters are monitored and documented automatically during the reprocessing cycle. A barcode scanner is used to record TEE probe identification and the user.

**Time efficiency of the reprocessing procedure**

With the exception of the Soluscope Serie TEE, the reprocessing procedure for the two other methods is purely manual (Tristel) or semi-automated (Optima C10). Manual reprocessing using the **Tristel Trio Wipes System** ties up staff resources as all reprocessing steps need to be performed and documented by the user. With the exception of the Tristel Sporidical Wipes, which are activated within 15 seconds by an activator foam, the impregnated wipes can be used immediately. A contact time of 30 seconds needs to be respected for disinfection.

With the **Optima C10**, a total time of 30 minutes is required for disinfection (20 minutes) and rinsing (10 minutes). Cleaning of the TEE probe needs to be performed separately prior to the automated reprocessing procedure. Depending on the specified product contact time in the immersion disinfection process, this can take up to 20 minutes with subsequent rinsing and drying, meaning a total time of up to 50 minutes is required for the entire reprocessing procedure.

<b>Tristel Trio Wipes System</b>	<b>Optima C10</b>	<b>Soluscope Serie TEE</b>
Visible contamination is wiped off using the Pre-Clean Wipe. Depending on the degree of contamination, cleaning may need to be repeated using multiple wipes. This step is documented in the log book.	The TEE probe is cleaned manually before it is loaded into the Optima C10. This involves preparing the probe with an enzymatic cleaner or disinfecting cleaning agent. It is then thoroughly rinsed with water of drinking water quality and wiped dry with a cloth. This step needs to be documented by the user.	The cleaning procedure is the first step of the automated process. This involves inserting the submersible distal part of the TEE probe into the cleaning/disinfection chamber and isolating it from the non-watertight handle during reprocessing by means of a retaining and sealing ring. Then the entire reprocessing procedure is started. Process monitoring and documentation are fully automatic.

**Tab. 1:** Cleaning.

Tristel Trio Wipes System	Optima C10	Soluscope Serie TEE
<p>Applying the activator foam generates chlorine dioxide on the Sporidical Wipe within 15 seconds. A chlorine-like odour confirms that the wipe is ready to use. After the probe has been wiped thoroughly, the required contact time of the disinfectant is 30 seconds. This step is documented in the log book.</p>	<p>The probe is disinfected in 20 minutes. This involves inserting the distal part of the probe into the immersion tube and connecting the TEE probe plug for the leak test. The disinfectant required for the process is automatically extracted from a media canister and then collected again in the same canister after disinfection. The reusable disinfection solution needs to be replaced by new solution at the end of the holding time or when the degree of contamination of the disinfectant (visible soiling) makes this necessary. An electrical residual current measurement (leak test) is performed during disinfection.</p>	<p>After the fully automatic cleaning process, the probe is gently disinfected for five minutes using a peracetic acid-based disinfectant developed by the manufacturer and adapted to the device. Documentation is fully automated. The electrical residual current measurement (leak test) is performed during the reprocessing procedure.</p>

Tab. 2: Disinfection.

Tristel Trio Wipes System	Optima C10	Soluscope Serie TEE
<p>Residual foam and disinfectant is wiped off using the Rinse Wipe. The wiping process should only be carried out once to prevent the wipe from fraying. This step is documented in the log book.</p>	<p>The fully automatic rinsing process takes place after disinfection and takes ten minutes. This is followed by automatic documentation of the 30-minute semi-automated reprocessing procedure via a printout.</p>	<p>The rinsing process takes place after disinfection and takes a few minutes. Documentation is fully automated.</p>

Tab. 3: Rinsing.

Tristel Trio Wipes System	Optima C10	Soluscope Serie TEE
<p>Instructions on drying TEE probes after rinsing are not provided.</p>	<p>Instructions on drying TEE probes after rinsing are not provided.</p>	<p>Drying takes place inside the reprocessing device. This is followed by automated documentation of the 14-minute reprocessing cycle.</p>

Tab. 4: Drying.

The entire reprocessing cycle including cleaning, disinfection, rinsing and drying takes a total of just 14 minutes with the **Soluscope Serie TEE**.

### Handling during TEE probe reprocessing procedures

#### **Tristel Trio Wipes (manual reprocessing)**

Manual reprocessing using the Tristel Trio Wipes System is based on impregnated wipes and solely involves wipe disinfection. Every step including cleaning, disinfection and rinsing needs to be meticulously executed and documented.

#### **Optima C10 (semi-automated reprocessing by immersion)**

TEE probe cleaning needs to be performed manually before the probes are loaded into the Optima C10. Hooking the distal part of the probe into the reprocessing tube and connecting the probe plug for the leak test is quick and easy. The TEE probe is disinfected and rinsed automatically by adding and draining the disinfectant and the reprocessed water for rinsing. The reusable disinfectant required for the process is extracted from a media canister and returned to the same canister after disinfection. The operating instructions recommend checking the minimum shelf life of the liquid disinfectant, its degree of contamination and the minimum quantity of 2 litres before each reprocessing operation. The standing time of the disinfectant specified by the manufacturer can be entered in the monitoring process parameters. Documentation of the semi-automated reprocessing procedure including the leak test is fully automatic.

#### **Soluscope Serie TEE (automated reprocessing)**

Cleaning, the leak test, disinfection, rinsing and drying steps are all performed fully automatically in just 14 minutes with the Soluscope Serie TEE. To load the device, the submersible distal part of the TEE probe is inserted into the cleaning/disinfection chamber and isolated from the non-watertight handle during reprocessing by means of a retaining and sealing ring. The automatic dispensing system removes the required quantity of process chemicals stored in the device for

the reprocessing procedure and disposes these via the water drain once reprocessing is complete. Process documentation including process monitoring is fully automatic.

### Protection of employees during disinfection

#### **Tristel Trio Wipes**

The Tristel Sporidical Wipe is as a virucidal oxidizing product based on chlorine dioxide under CLP classification not classified as a hazardous substance.

#### **Optima C10**

During the reprocessing procedure of the Optima C10, the use of aldehyde-based disinfectants can cause chemical vapours to evaporate, resulting in irritation of mucous membranes and allergic reactions. To prevent or minimize risks, the manufacturer's operating instructions recommend connecting the exhaust air from the reprocessing device to the room ventilation system or routing the air directly to the outside air area. Preparation of the ready-to-use disinfectant solution requires the general protective measures for handling disinfectant concentrates.

#### **Soluscope Serie TEE**

The use of peracetic acid as disinfectant and a pre-disinfecting cleaning solution during reprocessing completely eliminates the need for aldehyde-based process chemicals. The HEPA (High Efficiency Particulate Air) filtered exhaust air is free from odours and vapours. Please observe the instructions for safe handling of disinfectant concentrates when using peracetic acid.

### TEE probe transport and storage system

As part of the entire reprocessing chain, safe storage and protected transport in contamination-proof and sealed containers is also important to prevent recontamination of the TEE probe through improper storage or transport. In the event of recontamination, the reprocessing procedure needs to be repeated, which increases the time and costs involved for reprocessing a probe. Out of the three methods compared, Soluscope

is the only system that has a coordinated storage and transport solution<sup>7</sup> from the same provider. The TEE probes stored in the DSC 8000<sup>8</sup> drying and storage cabinet are recorded and documented using a barcode scanner via a touch screen display. This also facilitates traceability, documented storage time and resource planning for the stored TEE probes.

### Are the TEE probe reprocessing procedures for the three methods valid?

The German Society for Hospital Hygiene (DGKH) explicitly indicated in April 2017 that only validated methods may be used for reprocessing semi-critical medical devices. Assessing a process as valid requires proof of the constant efficacy of manual reprocessing methods. According to the judgement of the Higher Administrative Court (OVG) of North Rhine-Westphalia, ref. 13 A 24222/09 of 29.09.2010, this is not possible for critical medical devices due to the user's inability to ensure reproducibility of the individual reprocessing steps.

This concerns both manual reprocessing using the Tristel Trio Wipes System and semi-automated reprocessing using the Optima C10. Correct performance of reprocessing is not automatically guaranteed by these two methods as it is impossible to exclude human application errors. This means that the effectiveness of the method may not be reproducible and therefore cannot be considered valid. This also aligns with the Commission for Hospital Hygiene and Infection Prevention (KRINKO) recommendation of 2012, which stipulates that automated methods for cleaning and disinfection are easier to validate and that their use should be preferably used.

#### Tristel Trio Wipes

With Tristel Trio Wipes, there is a risk of incomplete cleaning and disinfection of the TEE probe, depending on the user's actions. For a full virucidal effect, the user needs to observe contact times when activating the wipe and during disinfection. Although documented standard work instructions and employee qualifications can be used to standardise process sequences, potential human application errors mean that it cannot be guaranteed that the procedure carried out is valid.


Fig. 4: Steps for wipe disinfection.

**Fig. 5:** Automated reprocessing using the Soluscope Serie TEE.

### **Optima C10**

Disinfection including a leak test and rinsing is automatic in the semi-automated reprocessing procedure of the Optima C10. Commercially available, authorised disinfectants can be used for TEE probe reprocessing. The manufacturer does not provide any optimised process chemicals for the reprocessing procedure. The user is responsible for using the correct concentration of disinfectant. The minimum shelf life and minimum quantity of the reusable disinfectant (2 litres) also need to be checked before starting the reprocessing procedure as per the user manual. If, for example, the disinfectant is heavily contaminated, it should be replaced before the end of its standing time. Manual pre-cleaning of medical devices and preparation and use of process chemicals need be performed according to documented standard work instructions. However, potential human application errors during pre-cleaning and use of the process chemicals mean that it cannot be guaranteed that the procedure carried out is valid.

### **Soluscope Serie TEE**

TEE probe reprocessing in the Soluscope Serie TEE is completely automated, including cleaning, disinfection, rinsing and drying. Every process step is documented automatically. During operation, there is constant monitoring of all process parameters to ensure adherence to nominal values. The Soluscope process chemicals are adapted to the reprocessing procedure, correctly dosed via the automatic dispensing system and disposed via the drain connection after use. Since all reprocessing steps are process-monitored, programme-controlled and fully automatic, the entire TEE probe reprocessing procedure of the Soluscope Serie TEE is valid.

The effectiveness of the cleaning result can therefore be automatically reproduced at any time. The integrated process monitoring and documentation of the entire reprocessing procedure ensures compliant, valid TEE probe reprocessing. The reprocessing procedure is easy to carry out, eliminates application errors and reduces staff costs during TEE probe reprocessing. The automatic leak test also helps to maintain the value of the probes as damage can be detected early.

## Summary

Reprocessing TEE probes is a major challenge in the daily routine of hospitals. Due to their high acquisition costs, TEE probes are often only available in small quantities, but they are also required for a high number of patients. This increases time pressure during reprocessing and can increase the risk of TEE probes being insufficiently cleaned and disinfected in the case of manual and semi-automated reprocessing. This is why it is particularly important that anyone entrusted with the task of reprocessing has the necessary expertise. The German Society for Hospital Hygiene (DKGH) explicitly indicated in April 2017 that only validated methods may be used for reprocessing semi-critical medical devices.

Assessing a process as valid requires proof of the constant efficacy of manual reprocessing methods. Even if the process is performed consistently in accordance with documented standard work instructions and suitable disinfectants adapted to the medical device are used, the potential for human error means that it cannot be guaranteed that the procedure carried out is valid. In the Soluscope Serie's fully automatic TEE probe reprocessing procedure, all reprocessing steps are process-monitored, programme-controlled and fully automatic. The effectiveness of the cleaning result can therefore be automatically reproduced at any time. The integrated process monitoring and documentation of the entire reprocessing procedure ensures compliant, valid TEE probe reprocessing. Even the process chemicals are provided in a ready-to-use state for use in the device. This makes it a safe procedure for the user and the patient.

## References

1. German Society for Sterile Supplies (DGSV): DGSV flow chart for classification of medical devices 2013; [https://www.dgsv-ev.de/wp-content/uploads/2016/09/AKQ\\_ZT\\_1\\_2013-3.pdf](https://www.dgsv-ev.de/wp-content/uploads/2016/09/AKQ_ZT_1_2013-3.pdf).
2. <http://www.krankenhaushygiene.de/informationen/624>.
3. Judgement of the Higher Administrative Court (OVG) of North Rhine-Westphalia, ref. 13 A 24222/09 of 29.09.2010; <https://openjur.de/u/144317.html>.
4. Soluscope video via Vimeo; <https://vimeo.com/152684309>.
5. Optima C10 brochure version 2011; [http://tee-reinigung.de/C10%20neue%20Version%20Prospekt\\_v2.pdf](http://tee-reinigung.de/C10%20neue%20Version%20Prospekt_v2.pdf).
6. Soluscope Serie TEE brochure; 07/2019; <https://www.soluscope.com/eng/page/washer-disinfector-for-probes-serie-tee-toe-0>.
7. Soluscope Transport brochure, version 03/2020; <https://www.soluscope.com/eng/page/storage-and-transport-solution>.
8. Soluscope DSC8000 brochure, version 03/2020; <https://www.soluscope.com/eng/page/endoscope-storage-drying-cabinet-dsc8000>.
9. RKI, Requirements for the reprocessing of medical devices. Recommendation issued by the Commission for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute (RKI) and the German Federal Institute for Drugs and Medical Devices (BfArM), German Federal Health Bulletin 2012.

# Reprocessing thermolabile endoscopes – the latest information and discussions

## Author

Dr Holger Biering  
 Holger Biering Consulting  
 Member of standardisation and  
 guideline groups for the reprocessing  
 of medical instruments in Germany,  
 Europe and the USA  
 41516 Grevenbroich, Germany  
 holger.biering@web.de

### Holger Biering

The improvement of workflows and processes for reprocessing thermolabile endoscopes with the aim of increasing patient safety during endoscopic examinations and procedures is a key topic featured in lots of technical literature in the

field and is the subject of working groups at national and international level with the intention of drawing up appropriate standards and recommendations. In the following, some of the results of this work will be presented and the current state of the discussion on this topic will be outlined.

### Sterilisation or disinfection

The issue of whether or not thermolabile endoscopes need to be sterilised is heatedly debated around the world. It is undisputed that low-germ (sterilised) endoscopes should be used in sterile areas of the body. There is also an ongoing discussion in Germany about the necessity to sterilise flexible bronchoscopes and cystoscopes.<sup>1,2,3</sup>

Based on a series of publications on infections with multi-resistant microorganisms associated with endoscopic retrograde cholangiopancreatography (ERCP) in various countries, the American standardisation organisation, the “Association for the Advancement of Medical Instrumentation” (AAMI), organised an expert forum where a general change to the risk classification of thermolabile endoscopes was discussed. According to this proposal, the classification should be changed as follows:

- Semi-critical medical devices: endoscopes which are using for diagnosis.

- Critical medical devices: endoscopes which are intended for therapeutic use (e.g. biopsies) or which may come into direct or indirect contact with non-intact skin or mucous membrane.<sup>4</sup>

According to this risk classification, the majority of endoscopes used in gastroenterology, for example, would have to be sterilised after each use on a patient. This change to the classification is justified by the increased safety margin of “overkill” sterilisation procedures compared to disinfection.<sup>5</sup> The aim was to include the results from the AAMI expert forum in corresponding American guidelines and standards.<sup>4</sup> In Europe and Germany, there are conflicting positions on this issue. Both the change to the risk classification and the transition from routine disinfection to sterilisation are not seen as conducive to increased patient safety for various reasons.<sup>6,7,8</sup>

In the USA, a revision of the corresponding AAMI Standard ST 91 “Flexible and semi-rigid endoscope processing in health care facilities”<sup>9</sup> was started after the expert forum and is currently nearing completion. A change to the risk classification is no longer included in the current version, but a gradual development from routine disinfection to sterilisation is recommended. It should be noted that in the USA, in addition to low-temperature sterilisation processes working with gaseous sterilisation media, such as ethylene oxide sterilisation and hydrogen peroxide sterilisation with and without plasma, it is also possible to use a liquid chemical sterilisation (LCS) process. LCS is a special type of disinfection based on the specific requirements for the approval of agents for final disinfection by the Food and Drug Administration (FDA) in the USA. The test requirements for agents for final disinfection at room temperature are different in Europe and the USA and not directly comparable. In simple terms, agents for final disinfection in the USA have to prove sporicidal efficacy in addition to the activity spectrum required in Europe.

When it comes to using the agents, a distinction is then made between:

- High level disinfection (HLD):  
Activity spectrum comparable with the requirements for final disinfection in Europe at a defined application concentration and for a defined application time.
- Liquid chemical sterilisation:  
Sporicidal efficacy in addition to the HLD activity spectrum at the same application concentration, but for a longer application time.

This means in the USA, any HLD disinfectant can also be used for LCS with an extended application time. With this extended application time, it can be assumed that the instrument materials are exposed to higher material stress and wear. Both disinfection processes are not relevant in Europe.

### Effectiveness of low-temperature sterilisation and disinfectants in the presence of inorganic and organic contamination

It is undisputed that only well-cleaned medical instruments can be successfully decontaminated in the subsequent disinfection and/or sterilisation reprocessing stages. This leads to the question of what cleaning level needs to be achieved. For a long time, a prescribed value of 6.4 µg/cm<sup>2</sup> was used as a guide value in the USA for the residual protein content on the instrument surface.<sup>10</sup>

Within the framework of the revision of Part 5 of the standard ISO 15883, criteria for tolerable residual content on instrument surfaces for protein as well as for other substances to be determined analytically, such as total organic carbon, carbohydrate, haemoglobin, ATP and endotoxin, are to be defined as a measure of purity.<sup>11</sup> The analytical criteria for protein are:

- Warning limit, ≥3 µg/cm<sup>2</sup>
- Intervention limit, ≥6.4 µg/cm<sup>2</sup>

In Europe, the effectiveness of instrument disinfectants is tested according to the standard BS EN 14885.<sup>12</sup> The tests for agents for final disinfection involve quantitative suspension tests and germ carrier tests with a low protein concentration using water of standardised hardness, i.e. inorganic and organic contamination are considered in the tests. At 16 µg/cm<sup>2</sup>, the applicable protein concentration in this process offers more than twice the safety margin of the intervention limit specified in the

standard ISO 15883-5. As part of the discussion on the sterilisation of bronchoscopes and cystoscopes, it was pointed out that publications from the late 1990s do not specify the effectiveness of low-temperature sterilisation procedures in the presence of inorganic and organic contamination.<sup>1</sup>

What is the situation 20 years later? To date, there is no international standard for hydrogen peroxide sterilisation that defines the special requirements for the process and the corresponding tests. A standard on this topic is currently being worked on at ISO level. The results of systematic studies into the effectiveness of an H<sub>2</sub>O<sub>2</sub> sterilisation process were presented at the Austrian Society for Sterile Supply (ÖGSV) congress in 2019. The test specimens used were made of metal and polyethylene non-woven fabric with and without organic contamination and were contaminated with spore suspension. The full effectiveness of the process could be confirmed for both carrier materials for the specimens without organic contamination. However, at a protein concentration of 20 µg/cm<sup>2</sup> and 10 µg/cm<sup>2</sup>, sufficient efficacy could not be determined for either of the carrier materials. It is the aim of further studies to determine whether efficacy can be granted with a protein concentration below the warning limit and intervention limit as well as with inorganic contamination.<sup>13</sup>

### Effectiveness of reprocessing against biofilm

The discussion on the effectiveness of reprocessing procedures for thermolabile endoscopes against biofilm is a key topic featured in lots of technical literature. In nature, biofilms are a complex mixture of different species that protect each other by forming intercellular substances. However, when testing the effectiveness of agents and procedures in laboratories, simplified single-species models are often used<sup>14,15</sup>, which are easier to inactivate or remove. It has been known for a long time that studies involving biofilms containing multiple species require a several times higher concentration of the disinfectant to inactivate the biofilm.<sup>16</sup> This statement was confirmed in a series of later studies by various authors. There are no known systematic studies on the effectiveness or testing of low-temperature sterilisation procedures against biofilm. This is not the case when it comes to machine-based cleaning/disin-

fection processes. When reprocessing in practice, it is a mandatory requirement that the removal of biofilm in the context of type testing of washer-disinfectors for reprocessing thermolabile endoscopes in accordance with ISO 15883-4 is tested by the manufacturer of the endoscope washer-disinfector.<sup>17,18</sup>

The method described in Annex B of this standard uses a multi-species biofilm model, which is more difficult to inactivate or remove than a single-species model. This means all endoscope washer-disinfectors for which a correct type test in accordance with the standard ISO 15883-4<sup>17</sup> was carried out, are also effective with regard to the elimination of biofilm. Additional measures to prevent the formation of biofilms and to eliminate them involve the effective implementation of manual cleaning steps, taking into account the current ESGE/ESGENA guidelines.<sup>6,7</sup> Mechanical support during cleaning or the use of enzymatic cleaners can be effective here.<sup>14,15</sup>

### Modified standard for washer-disinfectors

The standard ISO 15883-4 has been revised for washer-disinfectors for reprocessing thermolabile endoscopes.<sup>17</sup> Following a transitional period of 18 months to allow for the performance of extended tests within the scope of type testing by endoscope washer-disinfector manufacturers for devices already on the market, the standard has been in force since June 2020. Two major updates are of particular interest to users:

- Annex H describes a mandatory work instruction for endoscope washer-disinfector manufacturers within the scope of type testing for the endoscope washer-disinfector regarding the performance of tests for reprocessing individual endoscope types from different endoscope manufacturers. In order to keep the scope of testing at an acceptable level, endoscope manufacturers and endoscope washer-disinfector manufacturers can cooperate on this basis to group endoscopes into type test groups.
- Annex I provides an informative description of groups of endoscopes which can be used for validation in the context of performance qualification. This grouping of endoscopes is new for this standard, but corresponds to the “endoscope families” as per

Annex 3 of the German guideline on the validation of machine-based cleaning/disinfecting processes for reprocessing thermolabile endoscopes.<sup>19</sup>

In recent months, extensive work has been carried out by all endoscope washer-disinfector manufacturers to implement the requirements of Annex H. When purchasing a new endoscope washer-disinfector, the user can now ask detailed questions about whether the endoscopes used in their endoscopy unit are covered by the cleaning and disinfection tests within the scope of the type tests.

### Summary

Based on current information and discussions, the author proposes the following recommendations:

- Careful cleaning in accordance with the recommendations listed in the ESGE/ESGENA guidelines<sup>6,7</sup> is the prerequisite for successful reprocessing of thermolabile endoscopes in the subsequent process stages.
- Automatic cleaning and disinfection in an endoscope washer-disinfector is recommended for thermolabile endoscopes, which are classified as semi-critical medical devices, since this has been tested within the scope of type tests among others with regard to the effectiveness for removing biofilms and for a wide range of specific thermolabile endoscopes from various manufacturers.
- A disinfectant should be used in the endoscope washer-disinfector which has been tested in the presence of organic and inorganic contamination according to European test methods for application for the final disinfection. A high level of process reliability is achieved by the required scope of testing described above.
- Low-temperature sterilisation is only recommended for thermolabile endoscopes which are classified as critical medical devices in line with their medical application in accordance with the recommendation<sup>20</sup> by the Commission for Hospital Hygiene and Infection Prevention (KRINKO)/German Federal Institute for Drugs and Medical Devices (BfArM).

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Image p. 9: shutterstock

Circulation: 6.500

Publication schedule: three times a year

Printed on chlorine-free bleached paper

Only to be reprinted with the permission of the editorial team. Articles by named authors do not necessarily reflect the opinion of the editorial team. No liability is assumed for unsolicited manuscripts and photographs. The editorial team reserves the right to shorten letters from readers.

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Innovations Medical Vertriebs GmbH